

Small Area Sampling (SAS) Kit

The **Small Area Sampling (SAS) Kit** is designed to sample small areas using a Dacron swab or a gauze pad. Wet or dry samples can be taken with this kit. Two Dacron swabs are included in the SAS Kit that can be moistened with the buffer provided and used for wet areas or powders. The swab can be broken off into the buffer bottle, placed into the 50 mL centrifuge tube that is placed into a plastic bag for transport to a laboratory for analysis; or the SAS Kit can be used for taking dry samples using a gauze pad. The gauze is used to swipe the sample area and is then placed into the 50 mL centrifuge tube. The 50 mL tube is placed into a plastic bag for transport to a laboratory for analysis. A chain-of-custody form is included with each kit to record sample collection information. The SAS Kit can be used to test for munitions or munition fragments, terrorist lab/weapons material, packages with suspicious liquids or powders and recon of indoor areas. The SAS Kit is not intended to be used for testing soil, dirt or very dusty, dirty or porous surfaces. The kit is sterilized by Cobalt 60 irradiation at an ISO accredited facility. The sterility of each lot is also verified by random sterility testing.



Technical Specifications:

Size: 6" w x 11.5" h

Weight: 3 oz.

Buffer: 4mL of buffer

Kit Components:

- 1- Freestanding 50 cc tube
- 2- Dacron Swabs
- 1- 8 mL bottle, with dropper cap containing 4 mL buffer
- 1- Bitran®, liquid tight specimen bag
- 1- Chain-of-Custody form
- 1- Gauze pad
- 1- Instruction sheet

Advantages:

- Single use and easy-to-use
- Small, portable, and self-contained
- Suitable for sampling wet or dry surfaces
- Minimization of cross-contamination of samples
- Transportable following sample collection
- Sterile
- Hermetically sealed
- Traceable

Quantity discounts available-call for details

For Additional Information Contact:

TRADEWAYS LTD

184 Duke of Gloucester Street
Annapolis, Maryland 21401
United States of America

TEL: 1-410-295-0813

FAX: 1-410-295-0821

E-MAIL: Office@TradewaysUSA.com